

## **Rabbit anti-Sheep (Ovine) IgA Antibody (HRP)**

**Cat No:** LT6424742NIABU

**Antigen:** IgA

**Synonyms:** MB-1, IG-alpha, IGA, IgA, Igh-2

**Reactivity:** Sheep (Ovine)

**Host:** Rabbit

**Clonality:** Polyclonal

**Conjugate:** HRP

**Application:** ELISA

**Product Details Isotype:** IgG

**Characteristics:** Purified IgG conjugated to Horseradish Peroxidase (HRP)

**Application Notes:** Optimal working dilution should be determined by the investigator.

**Restrictions:** For Research Use only

**Handling Format:** Liquid

**Handling Concentration:** 1.0 mg/mL