

anti- 2019 nCOV N protein antibody

Product Information

Catalog No.:	LT80301MAB
Size:	100µg
Form:	liquid
Purification:	protein A purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	1G5
IsoType:	IgG1
Storage:	CB buffer, pH 7.5,-20°C for 12 months(Avoid repeated freeze / thaw cycles.)


Immunogen information

Immunogen:	2019 nCOV N protein
Observed MW:	48 kDa

Application

Reactivity:	2019 nCOV
Tested Application:	ELISA, WB, Colloidal Gold-Based
Recommended dilution:	ELISA: Use at an assay dependent concentration(Suggestions: Use a concentration of 0.1 - 1 µg/ml)

Image:



2019 nCOV N protein were subjected to SDS PAGE followed by western blot with anti- 2019 nCOV N protein Monoclonal antibody at dilution of 1ug/ml