

## **anti- 2019 nCOV S1 protein antibody**

### **Product Information**

Catalog No.:	LT82501MAB
Size:	100µg
Form:	liquid
Purification:	protein A+G purified
Purity:	≥95% as determined by SDS-PAGE
Host:	mouse
Clonality:	monoclonal
Clone ID:	2A6
IsoType:	IgG
Storage:	PBS pH 7.4,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)

### **Immunogen information**

Immunogen:	S1-RBD
Uniprot ID :	P0DTC2

### **Application**

Reactivity:	2019 nCOV S protein
Tested Application:	ELISA, WB
Recommended dilution:	ELISA: Use at an assay dependent concentration(Suggestions: Use a concentration of 0.1 - 1 µg/ml); WB: 1:1000-1:10000