

Anti-CD8 antibody (FITC)

Cat No: LT9052742NIABU

Overview Antigen: CD8

Synonyms: CD8, Leu2, MAL, p32

Reactivity: Cow (Bovine), Goat, Sheep (Ovine)

Host: Mouse **Clonality:** Monoclonal

Conjugate: FITC

Application: Flow Cytometry (FACS)

Immunogen: Ovine efferent lymphocytes Clone: 38-65 Isotype: IgG2a

Cross-Reactivity: Cow (Bovine), Goat

Characteristics: Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC)

Application Notes: Optimal working dilution should be determined by the investigator.

Restrictions: For Research Use only

Handling Format: Liquid

Concentration: 0.1 mg/ml