

MES, Potassium Salt

Product ID M5838

Description

MES potassium salt is a buffer that is most effective in the pH range of 5.5 to 6.7. It is the conjugate base of MES, free acid.

Product Information

Solubility	Water
Physical Form	Solid
Product Number	M5838
CAS Number	39946-25-3
Formula Weight	233.33
Formula	$C_6H_{12}KNO_4S$

Shipping and Storage

Storage Temp.	Room Temperature
Tariff Code	2934.99.9001

PhytoTechnology Laboratories®

14610 W 106th St. Lenexa, KS 66215

India Contact:

Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com