

Product Specification Sheet

Cat. #MA-30037

Mouse Monoclonal Anti Human Snail homolog 2 (SNAI2).; Ascites

SIZE : 100 ul

This gene encodes a member of the Snail family of C2H2-type zinc finger transcription factors. The encoded protein acts as a transcriptional repressor that binds to E-box motifs and is also likely to repress E-cadherin transcription in breast carcinoma. This protein is involved in epithelial-mesenchymal transitions and has antiapoptotic activity. The tumor suppressor protein p53 induces Slug expression in  $\gamma$ -irradiated cells; Slug protects damaged cells from apoptosis by repressing p53-induced transcription of the proapoptotic Bcl-2 family protein Puma. Mutations in this gene may be associated with sporadic cases of neural tube defects.

**Source of Antigen and Antibodies**

<b>Antigen</b>	Purified recombinant fragment of human SNAI2 expressed in E. Coli.
<b>Ab Host/type</b>	Balb/c mouse. IgG1 Ascetic fluid containing 0.05% sodium azide.
<b>2-Ab</b>	Goat Anti-mouse IgG-HRP conjugate Cat # 40320 (AP, biotin, FITC conjugates also available)
<b>-ve control IgG</b>	Cat # 20008-1, Mouse (non-immune) Serum IgG, purified, suitable for ELISA, Western, IHC as -ve control

**Isotype controls:**


Catalog#	Product Description
20102-101	Mouse IgG1 isotype control, purified
20102-101-APC	Mouse IgG1-APC conjugate (isotype)
20102-101-B	Mouse IgG1-Biotin conjugate (isotype)
20102-101-F	Mouse IgG1-FITC conjugate (isotype)
20102-101-FP	Mouse IgG1-FITC-PE conjugate (isotype)
20102-101-HP	Mouse IgG1-HRP conjugate (isotype)
20102-101-PC5	Mouse IgG1-PE-Cy5 conjugate (isotype)
20102-101-PE	Mouse IgG1-PE conjugate (isotype control)

**Suggested Dilutions:**

Western blot	1:500 – 1:2000
Immunohistochemistry (IHC):	N/A
Immunocytochemistry (ICC):	N/A
Flow cytometry (FCM):	N/A
ELISA	1:5000 – 1:100000

**Form:** Antibodies are supplied in PBS, pH 7.5, 0.05% azide and 0.1% BSA in liquid (0.5-1 mg/ml) or lyophilized in the same buffer. Reconstitute powder in 100 ul water or PBS. Store at -20°C or below is suitable size Aliquots.

**Shipping:** 4°C for solutions and room temp for powder.



**References:**

Biochem J. 2008 ;416(2):179-87. Mol Biol Cell. 2008 19(11):4875-87. Am J Pathol. 2009 174(6):2107-15.

*All products are for In vitro research use only.*

**Related material available from ADI**

MA-30037

150831V