

Product Specification Sheet

MOP4 (NPAS2) Antibodies

Cat. MOP41-S	Rabbit Anti-Human MOP4 antiserum # 1	SIZE: 100 ul
Cat. MOP41-A	Rabbit Anti-Human MOP4 IgG #1 (aff pure)	SIZE: 100 ug
Cat. MOP41-P	Human MOP4 Control/blocking peptide #1	SIZE: 100 ug

Several endogenous factors have been linked to rhythmicity or circadian behavior of living organisms. In *Drosophila*, the genes *period* (**dPer**) and *timeless* (*tim*), and in *Neurospora* *frequency* (*freq*), have been proposed to be responsible for their circadian rhythm. Recently human and mouse genes encoding a basic helix-loop-helix (bHLH) and Per-ARNT-Sim (PAS)-domain with significant similarity to the *Drosophila* Period have been reported. The cDNA sequences of hPER and mPer1 (also named *RIGU1*) are predicted to encode for proteins of length 1290 and 1291 amino acids respectively. Homologues of mPer1 designated **Per 2** (1257 aa) and **Per3** (1113 aa) have also been cloned. Both Per1 and Per2 levels show circadian rhythm in the SCN and eyes. It has been suggested that mPer regulates neuronal activity in the SCN. Mouse *Clock* (855 aa) is abundantly expressed in brain (SCN, pyramidal cortex, hippocampus) as well as in other tissues (eye, total brain, testis, ovaries, liver, heart, lung, and kidney). Although, *Clock* is constitutively expressed (not rhythmic) in the SCN, it may still be an important component of circadian machinery.

Basic-helix-loop-helix-PAS orphan protein, **MOP3** (Members Of PAS Superfamily; also known as **BMAL1/JAP3/PAS3**; 626 aa) is a general dimerization partner for several PAS superfamily of transcription regulators. MOP3 interacts with **MOP4** (also known as NPAS2; 626 aa), Clock, HIF1 α , and HIF2 α . MOP4 is a brain specific homolog of Clock.

Source of Antigen and Antibodies

Antigen	19-aa peptide of human MOP4 (protein accession #, refs 1) ; Designated (MOP41-P or control peptide) conjugated to KLH; epitope location ~ Mid-region
Ab Host/type	Rabbit, Polyclonal antiserum # MOP41-S and IgG, purified over antigen-agarose (Cat # MOP41-A). This antibody is also raised in chickens (cat # MOP42-S).
2-Ab	Goat Anti-rabbit IgG-HRP cat # 20320 (AP, biotin, FITC conjugates also available)
-ve control IgG	# 20009-1, Rabbit (non-immune) IgG, purified, suitable for ELISA, Western, IHC as -ve control

Form & Storage of Antibodies/Peptide Control

Antiserum (unpurified)	100ul solution	lyophilized powder
Supplied	0.05% azide, Reconstitute powder in 100 ul PBS	
Affinity pure IgG	100 ug/100ul solution	lyophilized powder

Supplied in **Buffer:** PBS+0.1% BSA
Reconstitute powder in PBS at 1 mg/ml

Control/blocking peptide
100 ug/100 ul solution lyophilized powder
Supplied in Buffer: PBS pH 7.5,
Reconstitute powder in PBS at 1 mg/ml.

Storage
Short-term: unopened, undiluted liquid vials at -20OC and powder at 4oC or -20oC..

Long-term: at -20C or below in suitable aliquots after reconstitution. Do not freeze and thaw and store working, diluted solutions.

Stability: 6-12 months at -20oC or below.
Shipping: 4oC for solutions and room temp for powder

Recommended Usage

Western Blotting (1:1K-5K for neat serum and 1-10 ug/ml for affinity pure using Chemiluminescence technique).

ELISA (1:10K-1:100K; using 50-100 ng of control peptide/well).

Histochemistry & Immunofluorescence: We recommend the use of affinity purified antibody at 2-20 ug/ml. This antibody has been cited in refs. 2.

Specificity & Cross-reactivity

Human MOP41-P sequence is 78% conserved in chicken and 73% in mouse MOP4/NPAS2. Antibody cross-reactivity with MOP4 from various species is not known. Control peptide, because of its low mol. Wt (<3 kDa), is not suitable for Western. It should be used for ELISA or antibody blocking experiments (use 5-10 ug control peptide per 1 ug of aff pure IgG or 1 ul antiserum) to confirm antibody specificity (see detailed protocol see detailed protocol at the web site).

General References:
Hogenesch JB et al (1997) J Biol. Chem. 272, 8581-8593;
Hogenesch JB et al (1998) PNAS 95, 5474; Zhou YD et al (1997) PNAS 94, 713-718.

Citations of ADI's MOP4 antibodies

Napoli C et al (2002) Circulation 105, 1360.

*This product is for *in vitro* research use only.

Related material available from ADI

Anti-Mouse/human Per1-3, Clock, MOP3-4; <i>Drosophila</i> Per, dClock, dBMAL, CRY1 and CRY2.	
MOP41-S-A-P	71210A