

Product Data Sheet

---

**Normal/Pre-immune Control Guinea Pig Serum, Mouse/G. pig tested**

**Cat#** NGPS-25-GM      **Size:**  25 ml      **Form:**  Liquid       Lyophilized

---

**Description**

A Donor population of adult Guinea Pig (either sex, mixed breed) is used to collect serum.

This lot of G. pig serum was also tested by ADI ELISA for species specific mouse proteins and found to be **negative** at 1:100 and 1:1000 serum dilution. Control mouse serum proteins generally gives highly positive ELISA reading at 1:100,000 and 1:1000,000 or more.

This lot of G. pig serum was also tested by ADI ELISA for species specific G. pig proteins and found to be **positive** at 1:10,000 and 1:100,000 serum dilution. Control G. pig serum generally gives highly positive ELISA reading at 1:100,000 and 1:1000,000 or more.

**Form and Storage**

Normal Guinea Pig serum is provided as liquid form. The product should be **stored frozen at -20°C** for short term storage. It is stable for a minimum of 1 year. Unless species otherwise, this preparation is not sterile.

**Suggested Uses**

This preparation of normal Guinea Pig serum is suitable for blocking in a given immunoassay or to serve as non-immune serum control for various primary antiserum. It may not be an ideal negative control for all primary antibodies as it is preferred to use the serum from the same Guinea Pig before any immunization. Control Guinea Pig serum should be diluted the same way as the primary antiserum. It may also be used to pre-treat the antibody beads during the immunoprecipitation to suppress non-specific adsorption of primary antibodies. Another use could be to use it as -ve control in ELISA for various primary antibodies, for which no pre-immune is available.

*This product is for in vitro research use only.*

**USDA Certification:** All animals are obtained from US sources and derived from USDA inspected and healthy stock.

**A Material Safety Data Sheet is not required for this product.**

The above product contain natural purified proteins in a buffer containing 0.02% sodium merthiolate or 0.05% azide. MSDS can be consulted for these two items. The product does not contain any hazardous components above 1% or any carcinogens above 0.1% as defined in 29 CFR 1910.1200, the OSHA Hazard Communication Standard. Please note that although the product is defined as not hazardous, it is still advisable to follow prudent laboratory practices when handling laboratory reagents. This material is sold for research purposes only and is not required to appear on the TSCA inventory. It is not intended for food, diagnostic, drug, household, agricultural or cosmetic use. Its use must be supervised by a technically qualified individual experienced in handling potentially hazardous chemicals.

The above information is correct to the best of our knowledge. Users should make independent decisions regarding completeness of the information based on all sources available. ADI shall not be held liable for any damage resulting from handling or contact with the above product.

---

Guinea Pig IgG, IgM, anti-Guinea Pig IgG and HRP, AP conjugates available.

Guinea Pig IgG, IgM ELISA kit available;

Guinea Pig IgG, IgM, IgA purified proteins

NGPS-25-GM

140320A

---

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)