

Product Data Sheet

**Normal/Pre-immune Control Pig/Swine Serum**

<b>Cat#</b>	NPIS-1	<b>Size:</b>	1 ml	<b>Form:</b>	Liquid	Lyophilized
<b>Cat#</b>	NPIS-10	<b>Size:</b>	10 ml	<b>Form:</b>	Liquid	Lyophilized
<b>Cat#</b>	NPIS-100	<b>Size:</b>	100 ml	<b>Form:</b>	Liquid	Lyophilized

**Description**

A healthy non-immunized population of Pig/Swine (Adult, mixed sex) was used to collect serum. Unless specified otherwise, this preparation is not sterile.

**Form and Storage**

The normal Pig/Swine serum is provided as liquid and 0.1% sodium azide as preservative or in lyophilized form. The product should be **stored at 4°C** for short term and **-20°C** for long term storage. It is stable for a minimum of 1 year.

The **lyophilized products** should be dissolved in PBS, pH 7.4 to bring it to original volume. Centrifuge the solution briefly (10,000 g, 15 min) to remove any particulate matter. It should be aliquoted and stored frozen for long term use.

**Suggested Uses**

This preparation of normal Pig/Swine serum is suitable for blocking in a given immunoassay or to serve as non-immune serum control for various primary antiserum. It may not be an ideal negative control for all primary antibodies as it is preferred to use the serum from the same Pig/Swine before any immunization. Control Pig/Swine serum should be diluted the same way as the primary antiserum. It may also be used to pre-treat the antibody beads during the immunoprecipitation to suppress non-specific adsorption of primary antibodies. Another use could be to use it as -ve control in ELISA for various primary antibodies, for which no pre-immune is available.

*This product is for in vitro research use only.*

**A Material Safety Data Sheet is not required for this product.**

The above product contain natural purified proteins in a buffer containing 0.02% sodium merthiolate or 0.05% azide. MSDS can be consulted for these two items. The product does not contain any hazardous components above 1% or any carcinogens above 0.1% as defined in 29 CFR 1910.1200, the OSHA Hazard Communication Standard. Please note that although the product is defined as not hazardous, it is still advisable to follow prudent laboratory practices when handling laboratory reagents. This material is sold for research purposes only and is not required to appear on the TSCA inventory. It is not intended for food, diagnostic, drug, household, agricultural or cosmetic use. Its use must be supervised by a technically qualified individual experienced in handling potentially hazardous chemicals.

The above information is correct to the best of our knowledge. Users should make independent decisions regarding completeness of the information based on all sources available. ADI shall not be held liable for any damage resulting from handling or contact with the above product.

Pig/Swine IgG, IgM, anti-Pig/Swine IgG and HRP, AP conjugates available.

Pig/Swine IgG, IgM ELISA kit available;

NPIS-1 -10 -100 -500

80709A

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi - 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)