

Product Specification Sheet

Mouse/Rat/Human OREXIN-A (Hypocretin-A) Fragment

Cat. OXA13-P Mouse/rat/human Orexin-A peptide 1-18 aa Fragment

SIZE:	100 ug	1 mg	5 mg	10 mg
FORM:	Soln	Lyophilized		

Several peptides associated with feeding behavior have been reported recently. Orexins (**Orexin-A** and **Orexin-B**) are a family of hypothalamic neuropeptides selectively expressed in the hypothalamus (1-2). Orexin-A and Orexin-B are derived from the same precursor (Prepro-orexin) by proteolytic cleavage. **Prepro-orexin** is 130 amino acid long peptide with a putative 33 AA secretory sequence, a hydrophobic core followed by residues with small polar side chains. The expression was detected in brain and to a small extent in testis (1-2). These neuropeptides bind and activate two closely related **Orexin receptors**--G-protein coupled receptors (GPCRs) **OX1R** and **OX2R**.

Source of peptide

A peptide sequence (**designated OXA13-P**) corresponding to 1-18 aa (<EPL PDC CRQ KTC SCR LYE, <E is Pyroglutamic acid; 2 disulfides at Cys6-Cys12 & Cys7-Cys14) of Orexin-A was synthesized and purified by HPLC. This sequence is the same in mouse, rat, and human Orexin-A.

Full length Orexin-A produced by us has been shown to be biologically active (3).

Form & Storage of Peptides

Control/blocking peptide

100 ug/100 ul solution	50 ug/50 ul lyophilized powder
Buffer: PBS, pH 7.5 and 0.05% sodium azide	
Reconstitute in the original vol. Of water	

Storage

Short-term: unopened, undiluted vials for less than a week at 4oC.

Long-term: at -20C or below in suitable aliquots after reconstitution. Do not freeze and thaw and store working, diluted solutions.

Stability: 6-12 months at -20oC or below.

Shipping: 4oC for solutions and room temp for powder.

General References:

Sakurai, T. et al. (1998) *Cell*, **92**, 573-585; DeLecca L et al (1998) *PNAS* **95**, 322-327; Kirchgeßner et al (1999) *Neuron* **24**, 941; Dube MG et al (1999) *Brian res.* **842**, 473-477; Jain MR et al (2000) *Regulatory Peptides* **87**, 18-24

*This product is for *in vitro* research use only.

Related material available from ADI:

Anti-Orexin 1 and 2 receptors, Anti-pre-proorexins

Anti-CART (Cocaine and Amphetamine related transcript)

Antibodies against Leptin, Leptin receptors, Tubby, Tub Homolog, Uncoupling proteins (**UCPs**) , **AGRP**, etc,

Melanocortin receptors, **BRS-3**, **Mahogany**

Western Blot recycling kit (Use the same blot to probe with multiple antibodies OX1, OX 2, etc.) **recycle blot at room temp in 5-10 min**; No mercaptoethanol or heating required).

Oxa13

rev. 40121S