

# Pipette Pump, Ergonomic with Thumbwheel, 2.0 mL

Product ID P302

## Description

Ergonomic pipette pump with a thumbwheel for easy dispensing. Compatible with pipettes up to 2.0 mL.

## Product Information

Capacity	2.0 mL
Use With	P990, P992
Product Number	P302

## Shipping and Storage

Storage Temp.	Room Temperature
Tariff Code	3926.90.9910

**PhytoTechnology Laboratories®**

14610 W 106th St. Lenexa, KS 66215

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)