

Pipette Pump, Ergonomic with Thumbwheel, 10.0 mL

Product ID P304

Description

Ergonomic pipette pump with a thumbwheel for easy dispensing. Compatible with 5.0 mL and 10.0 mL pipettes..

Product Information

Capacity	10.0 mL
Use With	P993, P994
Product Number	P304

Shipping and Storage

Storage Temp.	Room Temperature
Tariff Code	3926.90.9910

PhytoTechnology Laboratories®

14610 W 106th St. Lenexa, KS 66215

India Contact:

Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com