

Pipette Pump, Ergonomic with Thumbwheel, 25.0 mL

Product ID P306

Description

Ergonomic pipette pump with a thumbwheel for easy dispensing. Compatible with 25.0 mL pipettes

Product Information

Capacity	25.0 mL
Use With	P995
Product Number	P306

Shipping and Storage

Storage Temp.	Room Temperature
Tariff Code	3926.90.9910

PhytoTechnology Laboratories®

14610 W 106th St. Lenexa, KS 66215

India Contact:

Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com