

Picloram Solution (1 mg/mL)

Product ID P6717

Description

Picloram is an auxin that has been used as a 2,4-D substitute. This is the sterile filtered solution.

Product Information

Miscibility	Water
Physical Form	Liquid
Powder Form	P717
Sterility	Sterile Filtered
Product Number	P6717
CAS Number	1918-02-1
Formula Weight	241.46
Formula	$C_6H_3Cl_3N_2O_2$
Synonyms	4-Amino-3,5,6-trichloropicolinic Acid

Shipping and Storage

Storage Temp.	2 to 8 °C
Tariff Code	2942.00.0000

PhytoTechnology Laboratories®

14610 W 106th St. Lenexa, KS 66215

India Contact:

Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com