

Pectolyase Y-23

Product ID P686

Description

Pectolyase Y-23 is derived from *Aspergillus japonicas*. It is a purified maceration enzyme that contains two types of pectinase: endo-polygalacturonase and endo-pectin lyase. These two pectinases carry maceration stimulating factors that will macerate plant tissue when stimulated. Additionally, Pectolyase Y-23 when used in combination with cellulase, can isolate biologically active protoplasts.

Product Information

Solubility	Water
Physical Form	Powder
Product Number	P686
CAS Number	9033-35-6

Shipping and Storage

Storage Temp.	2 to 8 °C
Tariff Code	3507.90.7000

PhytoTechnology Laboratories®

14610 W 106th St. Lenexa, KS 66215

India Contact:

Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com