

Product Specification Sheet

Human Pancreatic Lipase (Triacylglycerol acylhydrolase)

Cat. #	PLIP15-N
Size	10 Units
Description	Human Pancreatic Lipase (partially pure) biologically active
Form	Solution in 0.1 M Tris containing 0.1 M NaCl and serine protease inhibitor.
concentration	~250 units/mg protein (see lot specific concn on the vial) One unit will liberate 1.0 μ mole of 2-monoglyceride from 1,2-diglyceride per min at pH 8.1 at 37 °C
Storage	Short term at 2-8oC; Long term at -20oC
Stability	~6 months year or more at 4oC; Avoid repeated freeze and thaw.

PLIP15-N

71221A