

Rabbit anti- *E. coli* O157:H7 (Polyclonal Antibody)

Target Protein:	<i>E. coli</i> O157:H7	
Catalog #:	410066	
Immunogen:	<i>E. coli</i> O157:H7	
Specificity:	reacts with <i>E. coli</i> O157:H7	
Clonality:	polyclonal	
Species:	Raised in rabbit	
Purity:	Protein-A Affinity purified	
Form:	0.1 mg (as determined by UV absorbance at 280nm) lyophilized, reconstitute in 1x PBS.	
Applications:	ELISA Dipsticks Immuno PCR	IMS Beads LAT Labeling
Storage:	Store at 2-8°C after dilution for up to 2 weeks. Aliquot and freeze at –20 °C to –80 °C for longer term storage. Avoid repeated freeze/thawing.	

For Laboratory research only, not for drug, diagnostics or other use.