

PathoScreen® Kit for Potato virus X (PVX)



Test Label: Alkaline Phosphatase
Test Format: Double Antibody Sandwich (DAS) ELISA
Capture Reagent: Polyclonal
Detection Reagent: Polyclonal

Intended Use:

This test can be used for detection of all strains of *Potato virus X (PVX)* in potato foliage, tuber sprouts and other solanaceous crops.

Includes:

- Antibody-coated 96-well microtiter plates
- Alkaline phosphatase enzyme conjugate
- RUB6, 1X Ready to use
- PNP substrate tablets, 5 mg each
- PNP substrate buffer, 5X concentrate
- Positive control
- PBST buffer, 20X concentrate, 50 mL
- Tween-20
- General extract buffer, GEB
- [User Guide](#)

PathoScreen® Kit, Potato virus X (PVX), 96 testwells	PSA 10002/0096
PathoScreen® Kit, Potato virus X (PVX), 288 testwells	PSA 10002/0288
PathoScreen® Kit, Potato virus X (PVX), 480 testwells	PSA 10002/0480