

PathoScreen[®] Kit for BPMV

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Intended Use:

The intended use of this test is to detect BPMV in plant parts, seed and beetles.

Includes:

- Antibody-coated 96-well microtiter plates
- Alkaline phosphatase enzyme conjugate
- ECI buffer, 5X concentrate
- PNP substrate tablets, 5 mg each
- PNP substrate buffer, 5X concentrate
- Positive control
- PBST buffer, 20X concentrate, 50 ml
- Tween-20
- General extract buffer, GEB
- [User Guide](#)



PathoScreen[®] Kit, Bean pod mottle virus (BPMV),PSA 46400/0096

PathoScreen[®] Kit, Bean pod mottle virus (BPMV),PSA 46400/0288

PathoScreen[®] Kit, Bean pod mottle virus (BPMV),PSA 46400/0480

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com