

## PathoScreen<sup>®</sup> Kit for AV-2

---

**Test Label:** Alkaline Phosphatase

**Test Format:** Antigen coated plate (ACP) ELISA

**Capture Reagent:** Polystyrene Plate

**Detection Reagent:** Polyclonal

### Intended Use:

---

This test is used to detect AV2 in asparagus.

### Includes:

- Detection antibody
- Alkaline phosphatase enzyme conjugate
- Indirect sample extraction buffer, 10X concentrate
- ECI buffer, 5X concentrate
- PNP substrate tablets, 5 mg each
- PNP substrate buffer, 5X concentrate
- Positive control
- 96 well microtiter plate, strip
- PBST buffer, 20X concentrate, 50ml
- [User Guide](#)



PathoScreen<sup>®</sup> Kit, Asparagus virus 2, PSA 71000/0096

PathoScreen<sup>®</sup> Kit, Asparagus virus 2, PSA 71000/0288

PathoScreen<sup>®</sup> Kit, Asparagus virus 2, PSA 71000/0480

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)