

PathoScreen[®] Kit for ZYMV

Test Label: Alkaline Phosphatase

Test Format: Double Antibody Sandwich (DAS) ELISA

Capture Reagent: Polyclonal

Detection Reagent: Polyclonal

Intended Use:

This is a qualitative serological assay for detection of *Zucchini yellow mosaic virus*(ZYMV) in plant foliage.

The test has been used for crop screens by Agdia Testing Services Department using cucurbit hosts, including squash, cucumber and melon.

The ZYMV ELISA detects the most important isolates from all over the world, including CT, USDA, SJBCA, CA, IT, NY, FL, Z18

Includes:

- Antibody-coated 96-well microtiter plates
- Alkaline phosphatase enzyme conjugate
- ECI buffer, 5X concentrate
- PNP substrate tablets, 5 mg each
- PNP substrate buffer, 5X concentrate
- Positive control
- PBST buffer, 20X concentrate, 50 ml
- Tween-20
- General extract buffer, GEB
- [User Guide](#)



PathoScreen[®] Kit, Zucchini yellow mosaic virus (ZYMV),PSA 77700/0096

PathoScreen[®] Kit, Zucchini yellow mosaic virus (ZYMV),PSA 77700/0288

PathoScreen[®] Kit, Zucchini yellow mosaic virus (ZYMV),PSA 77700/0480

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com