



Via Sant'Anna 131-135  
61030 Cartoceto PU (IT)  
Telephone + 39 (0)721830605  
Fax +39 (0)721837154  
e-mail: [info@diatheva.com](mailto:info@diatheva.com)  
[www.diatheva.com](http://www.diatheva.com)

**India Contact:****Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura,  
Delhi – 110034 (INDIA).  
Ph: +91-11-42208000, 42208111, 42208222  
Mobile: +91-9810521400  
Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com)  
Web: [www.lifetechindia.com](http://www.lifetechindia.com)

**Recombinant HIV-1 CRF02\_AG Tat****REP0059 - REP0060****100µg - 50µg****Description**

Multiple HIV-1 subtypes and circulating recombinant forms (CRFs) are known to cocirculate in Africa. In West Africa, the high prevalence of CRF02\_AG and cocirculation of subtype A, CRF01\_AE, CRF06\_cpx, and other complex intersubtype recombinants has been well documented. Cocirculation of CRF02\_AG strain with other genetic subtypes may give rise to novel and complex HIV-1 recombinants.

**Product type**

Recombinant protein


**Peptide**Recombinant protein: HIV-1Tat CRF02\_AG - **Accession number: AY371128** (truncated form)**Expression System***Escherichia coli***Tested by**

Western Blotting, SDS Page.

**Purity**

&gt; 90% pure estimated by SDS-PAGE (EU Ph. 5.0 § 2.5.31)

Purified recombinant CRF02\_AG tat (lane 1, molecular weight standard; lane 2, 1µg; lane 3, 2.5µg) was separated by SDS-PAGE (14% polyacrylamide) and stained with Coomassie Blue.

**Form**

Lyophilized.  
The protein should be reconstituted in apyrogenic sterile water or 1X PBS.

**Storage buffer**

Preservative 0.1% glycerol.

**Storage instructions**

Shipped at room temperature. The lyophilized protein is stable for 24 months if stored at -20°C. The reconstituted solution must be used immediately.