

Product Data Sheet

Cat # RP-1421

SARS-Associated Coronavirus Nucleocapsid Core Recombinant

Size: 100ug

SARS Coronavirus is an enveloped virus containing three outer structural proteins, namely the membrane (M), envelope (E), and spike (S) proteins. Spike (S)-glycoprotein of the virus interacts with a cellular receptor and mediates membrane fusion to allow viral entry into susceptible target cells. Accordingly, S-protein plays an important role in virus infection cycle and is the primary target of neutralizing antibodies. Immunoreactive with sera of SARS-infected individuals.

Source: The E.Coli derived 32 kDa recombinant protein contains the Nucleocapsid core protein 1-49 amino acids immunodominant regions. Purified by proprietary chromatographic technique

Applications and Suggested Dilutions: Protein is >95% pure as determined by 10% PAGE (coomassie staining). 50mM Tris-HCl, 60mM NaCl and 50% glycerol. Antigen in ELISA and Western blots, excellent antigen for detection of SARS with minimal specificity problems. Users must optimize the appropriate concentration and conditions for each assay.

Storage and Stability: Protein is shipped at ambient temperature. Upon arrival, store at -20°C. Five year frozen, 6 month at +4°C. If supplied in powder then reconstitute it in 100 ul water for 1 mg/ml stock and store in liquid at 4°C for ~1 week or aliquots in suitable size and store at -20°C for long term storage.

Usage:

This item is for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals

RP-1421 110131V

Alpha Diagnostic Intl Inc., 6203 Woodlake Center Dr, S an Antonio, T X 7 8 24 4 , U S A;

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400 Fax: +91-11-42208444 Email: customerservice@lifetechindia.com Website: www.lifetechindia.com