
Product Data Sheet

Cat # RP-1521

Influenza A Virus (H1N1) New Caledonia 20/99 IV 116

Size: 10 ug

Introduction:

H1N1 is a subtype specie of Influenza A virus. H1N1 Influenza Virus has mutated into various strains such as the Spanish Flu strain, mild human flu strains, endemic pig strains, and various strains found in birds. The Influenza A Virus is a globular particle about 100nm in diameter, sheathed in a lipid bilayer derived from the plasma membrane of its host. Studded in the lipid bilayer are two integral membrane proteins some 500 molecules of hemagglutinin ("H") and some 100 molecules of neuraminidase ("N"). Within the lipid bilayer are 3000 molecules of matrix protein and 8 pieces of RNA. Each of the 8 RNA molecules is associated with many copies of a nucleoprotein, several molecules of the three subunits of its RNA polymerase some "non-structural" protein molecules of uncertain function.

Description:

Allantoic fluid of 10 days old embryonated eggs, inoculated with influenza A virus, strain A/ New Caledonia/20/99 IVR 116. The Influenza Virus was purified by Ultracentrifugation with 10-40 % sucrose gradient. Thimerosal and beta propiolactone treatment. This product has been treated in a manner consistent with methods of inactivation. Generally accepted good laboratory practices appropriate to microbiological/viral safe handling practices and techniques are required at work.

Properties:

Greater than 90.0% as determined by:

Analysis by SDS-PAGE

Tested with anti-influenza A monoclonal antibodies in ELISA. Serological studies of influenza A virus, immunogen for antibody production.

Form

Sterile Filtered colorless solution. The H1N1 A/New Caledonia/20/99 IVR solution contains 0.1M NaCl, 10mM Tris Hcl, 1mM EDTA pH-8, 0.1 % sodium azide (NaN₃) and 0.005 % thimerosal.

Stability:

A/New Caledonia/20/99 IVR although stable 14 °C for 4 weeks, should be stored desiccated below -18 °C. **Please prevent freeze-thaw cycles.**

The product is for in vitro research use only.

RP-1521

90722A

Alpha Diagnostic Intl Inc., 6203 Woodlake Center Dr, San Antonio, TX 78244, U S A;

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400 Fax: +91-11-42208444 Email: customerservice@lifetechindia.com Website: www.lifetechindia.com