


Product Data Sheet

□ Cat # RP-1605 Recombinant (E.Coli) Dengue Virus Type 1 envelop protein (D-III), pure (>95%) **Size:** □ 50 ug

Dengue fever, also known as break bone fever, is an infectious tropical disease caused by the dengue virus. The term dengue fever came into general use only after 1828. The word dengue is Spanish for "affectation," "careful," or "fastidious. Dengue symptoms include fever, headache, muscle and joint pains, and a characteristic skin rash that is similar to measles. In a small proportion of cases the disease develops into the life-threatening dengue hemorrhagic fever (DHF), resulting in bleeding, low levels of blood platelets and blood plasma leakage, or into dengue shock syndrome, where dangerously low blood pressure occurs. Dengue is transmitted by several species of mosquito within the genus *Aedes*, principally *A. aegypti*. The virus has four different but related types 1-4 (Den1-4); infection with one type usually gives lifelong immunity to that type, but only short-term immunity to the others. Subsequent infection with a different type increases the risk of severe complications. There are up to 100 million cases of dengue fever worldwide every year; the most common occurrences are in urban parts of subtropical and tropical areas, such as Central and South America, parts of Africa, parts of Asia, the Caribbean and the Pacific.



Dengue fever virus (DENV) is an RNA virus of the family *Flaviviridae*. Other members of the same genus include yellow fever virus, West Nile virus, St. Louis encephalitis virus,

Japanese encephalitis virus, tick-borne encephalitis virus, Kyasanur forest disease virus, and Omsk hemorrhagic fever virus. Most are transmitted by arthropods (mosquitoes or ticks), and are therefore also referred to as arboviruses (arthropod-borne viruses). The dengue virus genome code for the three different types of structural proteins (**C**, **prM** and **E**) that form the virus particle and seven other types of non-structural proteins (**NS1**, **NS2a**, **NS2b**, **NS3**, **NS4a**, **NS4b**, **NS5**) that are only found in infected host cells and are required for replication of the virus. Dengue lab tests include cell culture, PCR, and antibody detection by ELISA. As there is no approved Dengue vaccine, prevention is sought by reducing the habitat and the number of mosquitoes and limiting exposure to bites.

FORMS & SOURCE:

Dengue Virus Type 1 envelop protein C-terminal domain 3 (D-III) was expressed as his-tag in *E. coli* and purified (>95%, ~15 kda). It is supplied in PBS, pH 7.4 and 0.05% azide (typically 0.5 mg/ml or lot specific concentration on the vial).

STORAGE & STABILITY:

Protein is shipped at ambient temperature. Upon arrival, Store at -20°C. 1 year frozen. One month in solution at room temperature. If supplied in powder then reconstitute it in 100ul water for 1mg/mL stock and store in liquid at 4°C for ~ 1week or aliquots in suitable size and store at -20°C for long term storage.

APPLICATION AND SUGGESTED DILUTIONS:

The expressed dengue envelope Type-1 domain III protein contains neutralizing epitopes and receptor binding domain.

References: Rodenhuis-Zybert IA (2010) *Cell Mol. Life. Sci.* 67, 2773-2786; Wiwanitkit, V (2010) Expert review of anti-infective therapy 8 (7): 841-5; Gubler DJ (1998) *Clin. Microbiol. Rev.* 11, 480-496; Villar L (2015) *NEJM* 372, 113-123

This item is for in vitro research use only (RUO).

Related Items

Catalog#	ProdDescription
540-100-DHG	Human Anti-Dengue Virus IgG ELISA kits, 96 tests
540-110-DHM	Human Anti-Dengue Virus IgM ELISA kits, 96 tests
540-120-DHG	Mouse Anti-Dengue Virus IgG ELISA kit, 96 tests
540-130-DHM	Mouse Anti-Dengue Virus IgM ELISA kit, 96 tests
540-RTH-100	Human Anti-Dengue Virus (NS1) rapid test, 100/pk
AB-14310	Mouse Anti-Dengue Type 2 (envelop) IgG, aff pure
AB-21120	Anti-Dengue Type 1-4 viruses antiserum
AB-21121	Monoclonal Anti-Dengue Virus Type 1-4 (pan, E antigen) IgG, culture medium
AB-21122	Monoclonal Anti-Dengue Virus Type 1-4 (pan, NS1 glycoprotein), culture medium
AB-21123	Monoclonal Anti-Dengue Virus Type 2, NS1 glycoprotein
AB-22122	# change to AB-21122; Monoclonal Anti-Dengue Virus Type 1-4 (pan, NS1 glycoprotein), culture medium
RP-1594	Recombinant (<i>E. coli</i>) Dengue Virus Type 1 E Antigen (DENV E), antigen grade (>95% pure)
RP-1595	Recombinant (<i>E. coli</i>) Dengue Virus Type 2 E Antigen (DENV E), antigen grade (>95% pure)
RP-1596	Recombinant (<i>E. coli</i>) Dengue Virus Type 3 E Antigen (DENV E), antigen grade (>95% pure)
RP-1597	Recombinant (<i>E. coli</i>) Dengue Virus Type 4 E Antigen (DENV E), antigen grade (>95% pure)
RP-1598	Recombinant (<i>E. coli</i>) Dengue Virus Type 2 (and epitopes for type 1, 2, and 3) E Antigen (DENV E) (23 kda, >95% pure)
RP-1599	Recombinant (<i>E. coli</i> , his-tag) Dengue Virus NS1, Type 4
RP-1600	Recombinant (<i>E. coli</i> , his-tag) Dengue Virus NS1, Type 3
RP-1601	Recombinant (<i>E. coli</i> , his-tag) Dengue Virus NS1, Type 1 RP-1602
RP-1602	Recombinant (<i>E. coli</i>) Dengue Virus Type 1-4, 2, and 3 envelop proteins antigen grade (>95% pure)
RP-1603	Recombinant (<i>E. coli</i>) Dengue Virus Type 3 envelop protein (D-III), pure (>95%)
RP-1604	Recombinant (<i>E. coli</i>) Dengue Virus Type 4 envelop protein (D-III), pure (>95%)
RP-1605	Recombinant (<i>E. coli</i>) Dengue Virus Type 1 envelop protein (D-III), pure (>95%)
RP-1606	Recombinant (<i>E. coli</i>) Dengue Virus Type 1-4 envelop+NS domains, pure (>95%)
RP-1607	Recombinant (<i>E. coli</i> , his-tag) Dengue Virus NS1 Type 2 immunodominant protein
RP-1608	Recombinant (<i>E. coli</i>) Dengue Virus Type 1 N-terminus envelop immunodominant regions, pure (>95%)
RP-1620	Recombinant (<i>E. coli</i> , his-tag) Dengue Virus NS1, Type 2 protein
RP-344	Recombinant Dengue Virus NS3 Type 1 protein
RP-345	Recombinant (<i>E. coli</i>) Dengue Virus NS1 c-end Type 2 protein
RP-346	Recombinant (<i>E. coli</i> , GST-tag) Dengue Virus NS1 n-end Type 2 protein
SP-100796-1	2A/2B Dengue Protease
SP-100797-1	2B/3, Dengue Protease Substrate
SP-100800-1	3/4A, Dengue Protease Substrate
RP-1601-Dengue-DENV-1-NS1-Protein.doc	

RP-1605-Dengue-1-Envelop-Protein 150819C

Alpha Diagnostic Intl Inc., 6203 Woodlake Center Dr, S an Antonio, T X 7 8 24 4 , U S A;

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400 Fax: +91-11-42208444 Email: customerservice@lifetechindia.com Website: www.lifetechindia.com