

Product Data Sheet

□ Cat # RP-388

Recombinant Human Matrix Metalloproteinase-8

Size: □ 5 ug

Full-length recombinant human neutrophil pro-collagenase (MMP-8), latent form. Matrix metalloproteinase 8 (MMP-8), or neutrophil collagenase, degrades interstitial collagens, acting preferentially on collagen type I. Increased full-length MMP-8 protein was associated with infiltration into the skin of neutrophils, which are the major cell type that expresses MMP-8. MMP-8 is synthesized and stored in specific granules in neutrophil leukocytes. MMP-8 activity is therefore regulated by factors such as surface-bound ligands (IgG or complement components) that release it through degranulation. Once released and activated through proteolytic or oxidative mechanisms, MMP-8 plays a major role in the connective tissue turnover that accompanies inflammatory processes.

**Source:** *Escherichia Coli*. Matrix Metalloproteinase-8 Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain having a molecular mass of 75 kDa. The MMP-8 is purified by proprietary chromatographic techniques. The protein Solution (100 µg/ml, 0.1 ml) in 0.05M Tris-HCl buffer, pH 7.5, containing 0.2M NaCl, 5mM CaCl<sub>2</sub>, 20µM ZnSO<sub>4</sub> and 0.1% BSA.

**Applications and Suggested Dilutions:** Greater than 90% as determined by SDS-PAGE. Used as a standard for analyzing mammalian collagenase activity. Users must optimize the appropriate concentration and conditions for each assay.

**Storage and Stability:** MMP-8 although stable at 14°C for 1 week, should be stored desiccated below -18°C. **Please prevent freeze-thaw cycles.** If supplied in powder then reconstitute it in 100 ul water for 1 mg/ml stock and store in liquid at 4°C for ~1 week or aliquots in suitable size and store at -20°C for long term storage.

**Biological Activity:** 100 units/ml after activation with APMA by solution assay method (10U=100 assays). One unit of collagenolytic activity is defined as the cleavage of 1 µg of collagen per minute.

**Usage:** This item is for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

RP-388

110208V

Alpha Diagnostic Intl Inc., 6203 Woodlake Center Dr, San Antonio, TX 78244, USA;

India Contact:

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400 Fax: +91-11-42208444 Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)