

Product Data Sheet

□ Cat # RP-404 Recombinant Dimethylglycine oxidase **Size:** □ 10 ug

Dimethylglycine oxidase (DMGO) is a covalent flavoenzyme from *Arthrobacter globiformis* that catalyzes the oxidative demethylation of dimethylglycine to yield sarcosine, formaldehyde, and hydrogen peroxide. The N-terminal region binds FAD covalently so it is yellowish.

Source: *Escherichia Coli*. Dimethylglycine oxidase Recombinant originated from *Arthrobacter globiformis* fused to His Tag at N-terminal produced in *E. Coli* is a single, non-glycosylated, polypeptide chain containing 850 amino acids and having a molecular mass of 92.1 kDa. The DMGO is purified by proprietary chromatographic techniques. The Recombinant Dimethylglycine Oxidase solution contains 20mM Tris-HCl pH7.5 and 20% glycerol.

Applications and Suggested Dilutions: Greater than 95.0% as determined by: (a) Analysis by RP-HPLC. (b) Analysis by SDS-PAGE. Users must optimize the appropriate concentration and conditions for each assay.

Storage and Stability: Dimethylglycine Oxidase Recombinant although stable at 4°C for 30 days, should be stored desiccated below -20°C for periods greater than 30 days. **Please prevent freeze-thaw cycles.** If supplied in powder then reconstitute it in 100 ul water for 1 mg/ml stock and store in liquid at 4°C for ~1 week or aliquots in suitable size and store at -20°C for long term storage.

Usage: This item is for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

RP-404 110208V

Alpha Diagnostic Intl Inc., 6203 Woodlake Center Dr, San Antonio, TX 78244, U S A;

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400 Fax: +91-11-42208444 Email: customerservice@lifetechindia.com Website: www.lifetechindia.com