
	Document ID:	TDS-RPL-001-100ML	Version:	003
	Date of Issue:	19-OCT-2020	Approved by:	Anna Becker
	Review Date:	10-JAN-2022	Signature:	
	Title:	<b>TECHNICAL DATA SHEET</b>		

### Papain Solution

**SPECIMEN**

Filtration	0.1µm sterile filtered
Product Code	RPL-001-100ML
Shelf Life	6 months from DOM
Storage Temperature	≤-15°C
Shipping Temperature	Dry ice

### Specifications

Physical and Chemical Analysis	Method	Specifications	Units
Appearance	Visual	Clear, colourless solution	n/a
pH @ 25°C	Electronic pH Meter	6.0 - 7.0	n/a
Osmolality	Osmometer	260 - 320	mOsm/kg
Endotoxin	LAL Kinetic	≤ 10.0	EU/ml
Sterility	Internally Validated	Pass	n/a

### Formulation

Components	Concentration (mg/L)
KCl	200.00
KH <sub>2</sub> PO <sub>4</sub>	200.00
NaCl	8000.00
Na <sub>2</sub> HPO <sub>4</sub>	1150.00
Papain	61.25

**Product Use:** Not intended for human or animal diagnostic or therapeutic use.

**THIS PRODUCT IS ANIMAL COMPONENT FREE**



India Contact:  
 Life Technologies (India) Pvt. Ltd.  
 306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
 Email: customerservice@lifetechindia.com Website: www.lifetechindia.com