

Product Data Sheet

□ Cat # RV26231-M Monoclonal Anti-M. tuberculosis Rv2031 dormant protein (H37Rv strain) IgG Size: 100 ul

Tuberculosis, MTB, or TB (short for tubercle bacillus) is a common, and in many cases lethal, infectious disease caused by various strains of mycobacteria, usually *Mycobacterium tuberculosis*. Tuberculosis typically attacks the lungs, but can also affect other parts of the body. Individuals with HIV are at risk for infection by tuberculosis due to their impaired immune system. The only currently available vaccine as of 2012 is bacillus Calmette–Guérin (**BCG** with live attenuated bacteria) which, while it is effective against disseminated disease in childhood, confers inconsistent protection against contracting pulmonary TB. A number of new TB vaccines are currently in phase I and II clinical trials. MVA85A (modified vaccinia Ankara 85A) is a subunit vaccine to BCG.

The closely related proteins of the **antigen 85 complex** are major secreted products of mycobacteria. Three closely related components, termed antigens **85A, 85B, and 85C**, have been demonstrated in *M. bovis* BCG and *M. tuberculosis*. Although the antigens are genetically distinct, they are highly homologous and cross-react with polyclonal and monoclonal antibodies raised against individual components. Sequence analysis revealed 85% identity between the *M. bovis* BCG 85A and 85B. Many other mycobacterial antigens have been identified, such as 71, 65, 38, 23, 19, 16, 14 and 12-kDa proteins. The **38-kDa** protein is an immunodominant lipoprotein antigen is specific only for the *M. tuberculosis* complex. It is the most extensively studied antigen. The 16-kDa antigen is an immunodominant antigen, frequently called 14 kDa, related to the family of low molecular weight heat-shock proteins. This antigen contains B-cell epitopes specific for the *M. tuberculosis* complex.

Antigen	Recombinant protein RV2623 of Dormancy regulon
Ab Host/type	Mouse, Monoclonal, IgG2a cat # RV26231-M) Supplied in PBS, pH 7.5 containing 0.05% azide,
2-ab	Goat Anti-mouse IgG-HRP conjugate Cat # 40320 (AP, biotin, FITC conjugates also available)
-ve control IgG	Cat # 20008-1, Mouse (non-immune) Serum IgG, purified, suitable for ELISA, Western, IHC as –ve control

Store antibody at 40oC for short term use and frozen at -20oC or below for long term storage. Reconstitute powder in water and store frozen.

Isotype Controls for mouse IgG2a

Catalog#	ProdDescription
20102-102	Mouse IgG2a isotype control, purified
20102-102-B	Mouse IgG2a-Biotin conjugate (isotype control)
20102-102-F	Mouse IgG2a-FITC conjugate (isotype control)
20102-102-FP	Mouse IgG2a-FITC-PE conjugate
20102-102-HP	Mouse IgG2a-HRP conjugate (isotype control)
20102-102-PC5	Mouse IgG2a-PE-Cy5 conjugate (isotype control)
20102-102-PE	Mouse IgG2a-PE conjugate (isotype control)

Specificity

RV26231-M Recognizes protein corresponding to RV2623 gene. Antibody crossreacts with *M. tuberculosis* RV3134 recombinant protein.

References: Cole ST (1998) Nature 393, 537-544; Renshaw PS (2005) EMBO J. 24, 2491-2498; Meher Ak (2006) FEBS J. 273, 1445-1462; Sorensen AL (1995) Inf. Immun. 63, 1710-1717; Berthet FX (1998) Microbiol. 144, 3195-3203; Brodin P (2006) Inf. Immun. 74, 88-98; Renshaw PS (2002) JBC 277, 21598-21603

This item is for LABORATORY RESEARCH USE ONLY.

Related Items

Catalog#	ProdDescription
HSP651-M	Monoclonal Anti-M. tuberculosis Heat Shock Protein 65 (hsp65/groEL-2/Cpn60-2) IgG
HSP651-P	Heat shock protein (<i>M. leprae</i> HSP65; 417-429) specific P62 peptide (LLQAAPALDKLKL, MW:1393.7)
HSP652-P	Heat shock protein (<i>M. leprae</i> / <i>M. tuberculosis</i> HSP65; 417-429) P38 peptide (AGGGVTLLQAAPALD, MW:1353.5)
HSP653-P	Heat shock protein (<i>M. leprae</i> HSP65; 343-355) P61 peptide (RVAQIRTEIENS, MW:1530.7)
HSP654-P	Heat shock protein (<i>M. bovis</i> HSP65; 243-255) indicator peptide in HLA-DQ2 binding assays (KPLLIAEDVEGEY, MW:1588.8)
HSP701-M	Monoclonal Anti-M. tuberculosis Heat Shock Protein 70 (hsp70/Dnak/ML2496) IgG
MTB061-M	Monoclonal Anti-Mycobacterium tuberculosis antigen (6kDa/ESAT-6) IgG
MTB06-R	Recombinant purified (<i>E. coli</i>) Mycobacterium tuberculosis antigen (6kDa/ESAT-6)
MTB161-M	Monoclonal Anti-Mycobacterium tuberculosis antigen (16kDa/HspX) IgG
MTB16-R	Recombinant purified (<i>E. coli</i>) Mycobacterium tuberculosis antigen (16kDa/HspX)
MTB381-M	Monoclonal Anti-Mycobacterium tuberculosis antigen 38kDa/Ag85B IgG
MTB38-R	Recombinant purified (<i>E. coli</i>) Mycobacterium tuberculosis antigen (38kDa/Ag85B)
MTB6381-S	Anti-M. Tuberculosis antigens (6Kda/ESAT+16kDa+38KDa/Ag85b mixture of proteins antiserum)
RP-628	Recombinant purified Mycobacterium Tuberculosis Heat Shock Protein 70 (hsp70/Dnak/ML2496)
RP-999	Recombinant purified Mycobacterium Tuberculosis major secretory protein Antigen 85B (38kda Antigen, Ag85b)
RV26231-M	130328A

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com