

Sterilization Indicator Strips, 3M Comply Plus, 8 inches

Product ID S7787

Description

Steam chemical indicator strips provide a visual if the autoclave's internal temperature reached sterilization (121 °C or 250 °F). Brand: 3M®

Product Information

Autoclavability	For use in steam autoclaves at 121 °C (250 °F)
Product Number	S7787
Length	8 inches (20.32 cm)

Shipping and Storage

Storage Temp.	Room Temperature
Tariff Code	3822.00.5090

PhytoTechnology Laboratories®

14610 W 106th St. Lenexa, KS 66215

India Contact:

Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com