

Schistosome Vaccine Related Proteins & Antibodies

(See Details at the website) http://4adi.com/commerce/catalog/spcategory.jsp?category_id=2818

| Items | Catalog # | Product Description | Product Type |
|---------------------|--|---|------------------|
| S. mansoni | SM231-R-50 | Recombinant (E. coli) Sm23 protein (S. mansoni, 1-218aa, His-tag >95%) | Rec. Protein |
| | SM251-R-100 | Sm25 antigenic peptide (S. mansoni, 36-82 aa, >95%, low endotoxin) | Rec. Protein |
| | SM281-R-50 | Recombinant (E. coli) Sm28/Smp28/GST28/GST-mu protein (S. mansoni, 1-211aa, His-tag >95%) | Rec. Protein |
| | SMAN11-M | Mouse monoclonal Anti-S. mansoni IgG for IF/ELISA | Mouse-mono |
| | SMP801-R-10 | Recombinant (E. coli) Sm-p80/Calpain/CANP (S. mansoni, 1-758aa, His-tag >95%) | Rec. Protein |
| | SMP801-S | Anti-Sm-p80/Calpain/CANP (S. mansoni, 1-758aa, His-tag >95%) antiserum | Rabbit-Poly |
| S. Japonicum | SM282-R-50 | Recombinant (E. coli) Sm28/Smp28/GST28/GST-mu protein (S. Japonicum, 1-211aa, no-tag >95%) | Rec. Protein |
| | SM282-S | Anti-Sm28/Smp28/GST28/GST-mu protein (S. Japonicum, 1-211aa) antiserum | Rabbit-Poly |
| | SM283-BTN | Monoclonal Anti-Sm28/Smp28/GST28/GST-mu protein (S. Japonicum) IgG-Biotin conj | Mouse-Mono |
| | SM283-M | Monoclonal Anti-Sm28/Smp28/GST28/GST-mu protein (S. Japonicum) IgG | Mouse-Mono |
| | SM283-R-50 | Recombinant (Sf9) Sm28/Smp28/GST28/GST-mu protein (S. Japonicum, 1-211aa, His-tag >95%, low endotoxin) | Rec. Protein |
| | GST11-A | Anti-Glutathione Transferase (GST, E. coli) IgG# 1, aff pure | Rabbit-Poly |
| | GST11-AP | Anti-Glutathione Transferase (GST, E. coli) IgG-AP conjugate | Rabbit-Poly |
| | GST11-C | Recombinant purified Glutathione Transferase, GST (E. coli), WB +ve control | Western control |
| | GST11-R | Purified Recombinant Glutathione Transferase (E. coli) protein | Rec. Protein |
| | GST11-S | Anti-Glutathione Transferase, GST (S. japonicum, E. coli) antiserum # 1 | Rabbit-Poly |
| | GST12-C | Recombinant purified Glutathione Transferase (GST)-Protein (S. japonicum, ~27 kda) control WB +ve control | Western control |
| | GST12-M | Monoclonal Anti-Glutathione Transferase, GST (S. japonicum, E. coli), ascites | Mouse-Mono |
| | GST13-A | Anti-Glutathione Transferase (GST, S. japonicum) IgG# 3 | Goat-Poly |
| | GST13-AS | Anti-Glutathione Transferase (GST, S. japonicum) IgG-Agarose (Aff matrix) | Goat-Poly |
| | GST13-BTN | Anti-Glutathione Transferase (GST, S. japonicum) IgG-Biotin Conjugate | Goat-Poly |
| | GST13-FITC | Anti-Glutathione Transferase (GST, S. japonicum) IgG-FITC Conjugate | Goat-Poly |
| | GST13-HRP | Anti-Glutathione Transferase (GST, S. japonicum) IgG-HRP Conjugate | Goat-Poly |
| | GST13-R | Purified Recombinant Glutathione Transferase-His(x6) tag (S. japonicum, GST-His) protein | Rec. Protein |
| | GST14-R | Purified Recombinant Glutathione Transferase-Ubiquitin fusion (S. japonicum, GST-Ub) protein | Rec. Protein |
| | GST15R-AS | Glutathione Transferase (GST, S. japonicum,) Protein-Agarose (Aff matrix) | Affinity support |
| GST16-BTN | Monoclonal Anti-Glutathione Transferase (GST, S. japonicum, E. Coli) IgG-Biotinylated | Mono-Mice | |
| GST16-M | Monoclonal Anti-Glutathione Transferase (GST, S. japonicum, E. Coli) IgG#2, purified | Mouse-Mono | |
| GST17-R | Purified Recombinant Glutathione Transferase (GST, S. japonicum) protein, native S. Japonicum (No tag) | Rec. Protein | |

Schistosoma-Vaccine-ELISA-Flr

160609A