
Cat# SP-100030-5

Description: [Tyr5,D-Trp6,8,9,Arg-NH210]-Neurokinin A (4-10)[Tyr5,D-Trp6,8,9,Arg-NH210]-Neurokinin A (4-10)

Sequence: Asp-Tyr-D-Trp-Val-D-Trp-D-Trp-Arg-NH₂; MW 1109.3

Size: 5 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-100030-5 101110A