
Cat# SP-100032-5

Description: [Pro7]-Neurokinin B[Pro7]-Neurokinin B

Sequence: Asp-Met-His-Asp-Phe-Phe-Pro-Gly-Leu-Met- NH₂; MW 1208.43

Size: 5 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-100032-5 101110A