
Cat# SP-100033-5

Description: [D-Pro2,D-Trp6,8,Nle10]-Neurokinin B[D-Pro2,D-Trp6,8,Nle10]-Neurokinin B

Sequence: Asp-D-Pro-His-Asp-Phe-D-Trp-Val-D-Trp-Leu-Nle-NH₂; MW: 1326.53

Size: 5 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-100033-5 101110A