

Cat# SP-100036-1

Description: Prepro-Neuromedin S (70-103) (human) (AA: Phe-Leu-Phe-His-Tyr-Ser-Arg-Thr-Gln-Glu-Ala-Thr-His-Pro-Val-Lys-Thr-Gly-Phe-Pro-Pro-Val-His-Pro-Leu-Met-His-Leu-Ala-Ala-Lys-Leu-Ala-Asn) (MW: 3857.5)

Size: 1 mg

Purity: >95%

Store: Desiccated at -20oC.

SP-100031-5 101110A