

Cat# SP-100040-1

Description: [Leu116]-Prepro-Neuromedin U (104-136) (human)[Leu116]-Prepro-Neuromedin U (104-136) (human)

Sequence: Phe-Leu-Phe-His-Tyr-Ser-Lys-Thr-Gln-Lys-Leu-Gly-Leu-Ser-Asn-Val-Val-Ser-Ser-Val-Val-His-Pro-Leu-Leu-Gln-Leu-Val-Pro-His-Leu-His-Glu; MW: 3768.45

Size: 1 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-100040-1 101110A