

Cat# SP-100047-1

Description: Adrenomedullin (11-50) (rat)Adrenomedullin (11-50) (rat)

Sequence: Ser-Thr-Gly-Cys-Arg-Phe-Gly-Thr-Cys-Thr-Met-Gln-Lys-Leu-Ala-His-Gln-Ile-Tyr-Gln-Phe-Thr-Asp-Lys-Asp-Lys-Asp-Gly-Met-Ala-Pro-Arg-Asn-Lys-Ile-Ser-Pro-Gln-Gly-Tyr-NH₂
(Disulfide bridge:Cys4-Cys9) (MW: 4521.2)

Size: 1 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-100047-1 101110A