

Cat# **SP-100050-1**

Description: Proadrenomedullin (45-92) (human) (AA: Glu-Leu-Arg-Met-Ser-Ser-Ser-Tyr-Pro-Thr-Gly-Leu-Ala-Asp-Val-Lys-Ala-Gly-Pro-Ala-Gln-Thr-Leu-Ile-Arg-Pro-Gln-Asp-Met-Lys-Gly-Ala-Ser-Arg-Ser-Pro-Glu-Asp-Ser-Ser-Pro-Asp-Ala-Ala-Arg-Ile-Arg-Val) (MW: 5114.76)

Size: 1 mg

Purity: >95%

Store: Desiccated at -20oC.

SP-100050-1 120703A