
Cat# SP-100062-1

Description: [Ala31, Aib32]-Neuropeptide Y (porcine)[Ala31, Aib32]-Neuropeptide Y (porcine)

Sequence: Tyr-Pro-Ser-Lys-Pro-Asp-Asn-Pro-Gly-Glu-Asp-Ala-Pro-Ala-Glu-Asp-Leu-Ala-Arg-Tyr-Tyr-Ser-Ala-Leu-Arg-His-Tyr-Ile-Asn-Leu-Ala-Aib-Arg-Gln-Arg-Tyr-NH₂; MW: 4194.63

Size: 1 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-100062-1 101110A