

**Cat#** SP-100063-1

**Description:** [Leu31,Pro34]-Neuropeptide Y (porcine)[Leu31,Pro34]-Neuropeptide Y (porcine)

**Sequence:** Tyr-Pro-Ser-Lys-Pro-Asp-Asn-Pro-Gly-Glu-Asp-Ala-Pro-Ala-Glu-Asp-Leu-Ala-Arg-Tyr-Tyr-Ser-Ala-Leu-Arg-His-Tyr-Ile-Asn-Leu-Leu-Thr-Arg-Pro-Arg-Tyr-NH<sub>2</sub>; MW: 4240.8

**Size:** 1 mg

**Purity:** >95%

**Store:** Desiccated at -20oC.

*All peptides are for in vitro research use only.*

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.  
[http://www.4adi.com/commerce/info/showpage.jsp?page\\_id=1088&category\\_id=2427](http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427)

SP-100063-1                    101110A