
Cat#	SP-100076-5
Description:	Ac-[Leu28,31]-Neuropeptide Y (24-36)Ac-[Leu28,31]-Neuropeptide Y (24-36)
Sequence:	Ac-Leu-Arg-His-Tyr-Leu-Asn-Leu-Leu-Thr-Arg-Gln-Arg-Tyr-NH ₂ (MW: 1787.1)
Size:	5 mg
Purity:	>95%
Store:	Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-100076-5 101110A