

Cat# SP-100086-5

Description: [Ser4,Ile8]-Oxytocin[Ser4,Ile8]-Oxytocin

Sequence: Cys-Tyr-Ile-Ser-Asn-Cys-Pro-Ile-Gly-NH2
(Disulfide bridge: Cys1-Cys6□; MW 966.15)

Size: 5 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-100086-5 101110A