
Cat# SP-100292-5

Description: [Cys(Acm)20,31]-EGF (20-31)[Cys(Acm)20,31]-EGF (20-31)

Sequence: Cys(Acm)-Met-His-Ile-Glu-Ser-Leu-Asp-Ser-Tyr-Thr-Cys(Acm); MW: 1543.77

Size: 5 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-100292-5 101110A