

Product Data Sheet

Cat# SP-100370-25

Description: [Sar1,Val5,Ala8]-Angiotensin II[Sar1,Val5,Ala8]-Angiotensin II

Sequence: Sar-Arg-Val-Tyr-Val-His-Pro-Ala; MW 912.1

Size: 25 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-100370-25 101110A