

Product Data Sheet

Cat# SP-100451-1

Description: [Leu31,Pro34]-Peptide YY (human)[Leu31,Pro34]-Peptide YY (human)

Sequence: Tyr-Pro-Ile-Lys-Pro-Glu-Ala-Pro-Gly-Glu-Asp-Ala-Ser-Pro-Glu-Glu-Leu-Asn-Arg-Tyr-Tyr-Ala-Ser-Leu-Arg-His-Tyr-Leu-Asn-Leu-Leu-Thr-Arg-Pro-Arg-Tyr-NH₂; MW: 4292.9

Size: 1 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-100451-1 101110A