

Cat# **SP-100456-1**

Description: PACAP-38 (16-38) (human, chicken, mouse, ovine, porcine, rat) (AA: Gln-Met-Ala-Val-Lys-Lys-Tyr-Leu-Ala-Ala-Val-Leu-Gly-Lys-Arg-Tyr-Lys-Gln-Arg-Val-Lys-Asn-Lys-NH₂) (MW: 2720.37)

Size: 1 mg

Purity: >95%

Store: Desiccated at -20oC.

SP-100456-1 120703A