

Product Data Sheet

Cat# **SP-100511-1**

Description: β -Defensin-3, human β -Defensin-3, human

Sequence: Gly-Ile-Ile-Asn-Thr-Leu-Gln-Lys-Tyr-Tyr-Cys-Arg-Val-Arg-Gly-Gly-Arg-Cys-Ala-Val-Leu-Ser-Cys-Leu-Pro-Lys-Glu-Glu-Gln-Ile-Gly-Lys-Cys-Ser-Thr-Arg-Gly-Arg-Lys-Cys-Cys-Arg-Arg-Lys-Lys (Disulfide bridge: Cys11-Cys40, Cys18-Cys33, Cys23-Cys41) (MW: 5155.22)

Size: 1 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-100511-1 101110A