
Cat# **SP-100515-5**

Description: [D-Arg8] - Dynorphin A (1 - 13), porcine[D-Arg8] - Dynorphin A (1 - 13), porcine

Sequence: Tyr-Gly-Gly-Phe-Leu-Arg-Arg-D-Arg-Arg-Pro-Lys-Leu-Lys; MW: 1647.01

Size: 5 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-100515-5 101110A