

Product Data Sheet

---

**Cat#**                **SP-100825-1**

**Description:**    Adrenomedullin (1- 52), porcineAdrenomedullin (1- 52), porcine

**Sequence:**        (AA: see antigen, sequence too long) (MW: 5971.8)

**Size:**                1 mg

**Purity:**             >95%

**Store:**              Desiccated at -20oC.

---

*All peptides are for in vitro research use only.*

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.  
[http://www.4adi.com/commerce/info/showpage.jsp?page\\_id=1088&category\\_id=2427](http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427)

**SP-100825-1**                101110A