

Cat# SP-100882-1

Description: [D-Ser14] - Humanin (HN)[D-Ser14] - Humanin (HN)

Sequence: Met-Ala-Pro-Arg-Gly-Phe-Ser-Cys-Leu-Leu-Leu-Leu-Thr-D-Ser-Glu-Ile-Asp-Leu-Pro-Val-Lys-arg-Arg-Ala; MW: 2687.28

Size: 1 mg

Purity: >95%

Store: Desiccated at -20oC.

All peptides are for in vitro research use only.

Please consult "Frequently asked questions" section at our website for Guidance on storage and solubility of the peptides.
http://www.4adi.com/commerce/info/showpage.jsp?page_id=1088&category_id=2427

SP-100882-1 101110A